

HOME-BASED NON-INVASIVE LUNG COMPLIANCE MONITORING DEVICE

Enable non-invasive measuring of lung compliance to monitor and characterize diseases that alter lungs' mechanical properties.

EDWARD CASTILLO

Associate Professor, Department of Biomedical Engineering, Cockrell School of Engineering, UT Austin

HSU-TING KUO

PhD Candidate, Department of Biomedical Engineering, Cockrell School of Engineering, UT Austin

- Over 16 million U.S. adults have lung diseases affecting the pulmonary parenchyma, costing over \$52 billion annually^{1,2,3}.
- Lung compliance (LC)** measures lung's ability to expand at a given pressure, serving as a key indicator of disorders that alter lung's mechanical properties.


	Common pulmonary parenchymal disorders		
	COPD (Chronic Obstructive Pulmonary Disease) ^{2, 4}	ILD (Interstitial Lung Disease) ^{3,5}	RP (Radiation pneumonitis) ⁶
Prevalence/ Incidence	16M (diagnosed) 24M (undiagnosed)	200K	15-40% of thoracic radiotherapy patients
Cost (USD annually)	50B	2B	10M
5-year survival rate (%)	40-70%	56% (Overall) 41% (IPF) ⁵	-
ER visits (in the US annually)	1.5M	13.7%-19.4% of cases	-

- Current lung compliance measuring tools are either **invasive** or **involve radiation exposure from CT imaging**.
- Clinical Unmet Need: No **non-invasive, home-based** solution for regular lung disease monitoring is currently available on the market.

Current LC measuring tools:	Non-invasive?	Precise, consistent pressure control	Frequent monitoring	One-time setup cost	Cost per exam
Mechanical ventilation ^{1, 2}		✓		\$2.5K – 5K	\$1,522
Pleural manometry ^{2, 3}				N/A	\$1,038-1,189 (Thoracentesis)
Esophageal manometry ⁴				N/A	\$924-1,391
4DCT imaging ^{5,6}	✓			\$300K -3M	\$6,773
Our Solution	✓	✓	✓	\$700*	\$ 0.75


* USD 700 covers the cost of the device and the mouthpiece

- Solution: A **non-invasive** system that calculates lung compliance (LC), an important metric in overall lung health, in an **easy-to-use, at-home** package that can be set up and used in minutes.
- Current developmental statue: **Prototyping**
 - Prototype **physical design**
 - Fine-tune **algorithms** for pressure regulation and measurement



- $LC = \frac{\Delta V}{\Delta P} = \frac{\text{Lung volume change}}{\text{Pressure change}}$ ← Calculated using our numerical pipeline
 ← $P_{control}$ across different time points
- Numerical pipeline to estimate airflow rate at all points in the breath cycle
- PID control to regulate pressure during normal breathing cycles

Numerical pipeline



Poiseuille's formula:


$$Q(t) = \frac{P_{control}(t) - P_{measure}(t)}{Res}$$

Numerical integration:

$$Q(t) \cong \frac{dV(t)}{dt} = \frac{V(t + \Delta t) - V(t - \Delta t)}{2\Delta t}$$


Lung compliance calculation:

$$LC = \frac{\Delta V}{\Delta P} = \frac{V(t_1) - V(t_2)}{P_{control}(t_1) - P_{control}(t_2)}$$



- Our team has demonstrated significant difference in 4DCT-derived lung compliance (LC) in **healthy** and **IPF** patients¹.
- We observed **significant decrease** in 4DCT-derived lung compliance over 6 months.

IPF patient



Non-IPF patient

